

# How Do I ... ? Or Mitosis Visualizer in brief

Makes high-res screenshot and stores it in an out folder

Starts or stops video capturing, video will be stored in an out folder

Use mouse to perform the following actions:  
L Button – select an object  
R Button – rotate camera  
R Button + Wheel - zoom

Rewinds simulation




Shows window with detailed information about selected object

Play

Suspends or resumes playback

Sets visualization speed as the model seconds to a real second ratio

Defines how objects will be visualized (colors, lighting models, etc)



Defines the name prefix from which image and video files will be started

Defines the directory in which images and videos will be stored

Defines the video storage format

Defines the image storage format

